

BATHLESS



DISINTEGRATION TESTER SYNCHRONOUS AND NON-SYNCHRONOUS (BATHLESS)

ULD-2SABL (RPS) AND ULD-2SABL (RPNS)

- USP, IP, EP Compliance.
- Synchronous drive with continuous timer and Non – Synchronous with Individual Timer.
- Smooth transition drive.
- Individual Temperature Probe.
- Individual Tablet Disintegration time registration.
- Bottom illumination for better visibility.
- Magnetic Basket for easy basket handling.
- Basket Raise park from beaker with Beep after test completion.
- Print Test Report & Validation Report.

KEY FEATURES OF BATHLESS DISINTEGRATION TESTER

TIME: Achieving temperature in Bathless Disintegration Tester is 50% faster than the Disintegration tester with Bath.

ELECTRICITY: As it achieves the temperature faster which helps in less consumption of Electricity.

USER-FRIENDLY: Easy for visual observation of Disintegration of tablet & for cleaning process also.

WATER: As it is Bathless , the usage of water become less.